

## Standard Operating Procedures: Preparing Bulk Samples for Radiocarbon Dating

Lab Name: Rosenheim Lab

Lab Location: Marine Science Laboratories, CMS, USF

Room #: 227A

Principal Investigator: Brad E. Rosenheim

Page: 1 of 2

Time Frame: 2 half days – can't be completed in a single day because samples must be combusted

<p><b>1. Process</b></p>	<ol style="list-style-type: none"><li>1. While wearing gloves, weigh out enough sediment to produce at least 24 <math>\mu\text{mol}</math> of <math>\text{CO}_2</math> and put sediment into a prepped quartz tube<ul style="list-style-type: none"><li>- See tube prep protocol</li></ul></li><li>2. Place 1-cm of quartz wool over the sediment to prevent the sediment from entering the vacuum line when the tube is evacuated</li><li>3. Attach a 1/2 - 1/4" fritted ultratorr to the inlet on the vacuum line and attach the quartz tube</li><li>4. Evacuate the tube of atmospheric pressure<ul style="list-style-type: none"><li>- Evacuation can take anywhere from 30 minutes to 6+ hours</li><li>- Check for leaks before flame sealing. Close the valve to the vacuum line and look for an increase in pressure. See leak testing information page.</li></ul></li><li>5. Once evacuated, flame seal the quartz tube and remove it from the line</li><li>6. Combust the tube at 900°C for four hours<ul style="list-style-type: none"><li>- See furnace protocol</li></ul></li><li>7. After the tube has cooled, attach the tube cracker to the inlet using a fritted ultratorr and crack open the tube</li><li>8. Move the <math>\text{CO}_2</math> to the next trap using liquid nitrogen</li><li>9. Move the sample down the vacuum line, alternating liquid nitrogen and liquid nitrogen/isopropanol slush to purify the sample of water vapor</li></ol>
--------------------------	---

	<p>10. Once the CO<sub>2</sub> reaches the Inficon, let the gas warm up to room temperature and record a pressure reading</p> <p>11. Get an additional pressure reading with slush</p> <p>12. Move the CO<sub>2</sub> back down the vacuum line, alternating liquid nitrogen and liquid nitrogen/isopropanol slush</p> <p>13. Transfer the gas into an empty, pre-combusted pyrex tube using liquid nitrogen and flame seal</p> <ul style="list-style-type: none"> <li>- See furnace protocol</li> </ul>
<b>2. Hazardous Chemical Class of Hazardous Chemical</b>	N/A
<b>3. Personal Protective Equipment</b>	Closed toed shoes, flame sealing glasses, hair pulled back
<b>4. Engineering/Ventilation Controls</b>	N/A
<b>5. Special Handling Procedures Storage Requirements</b>	N/A
<b>6. Spill Containment Accident Procedures</b>	N/A
<b>7. Waste disposal</b>	Broken tubes are disposed of in the cardboard glass disposal bin under the furnace
<b>8. Special Precautions Animal Use</b>	N/A
<b>9. Required Approvals</b>	N/A
<b>10. Decontamination</b>	Pre-combustion of quartz and glassware
<b>11. Designated Areas</b>	N/A